

## Product datasheet for **TA386833**


### CALML5 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Calmodulin-like protein 5 protein (1-146AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9NZT1</a>
Background:	Binds calcium. May be involved in terminal differentiation of keratinocytes.



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded human skin tissue using TA386833 at dilution of 1:100



Western blot  
All lanes: Calmodulin-like protein 5 antibody at 12µg/ml + A431 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 16 kDa  
Observed band size: 16 kDa