

## Product datasheet for **TA386905**

### **DCTN6 Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC, IP, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IP:1:200-1:2000
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Dynactin subunit 6 protein (1-190AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">O00399</a>
Background:	centrosome, cytosol, antigen processing and presentation of exogenous peptide antigen via MHC class II, ER to Golgi vesicle-mediated transport




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human gastric cancer using TA386905 at dilution of 1:100




Immunoprecipitating DCTN6 in A549 whole cell lysate


Lane 1: Rabbit control IgG instead of (1 $\mu$ g) instead of TA386905 in A549 whole cell lysate. For western blotting, a HRP-conjugated light chain specific antibody was used as the Secondary antibody (1/50000)

Lane 2: TA386905 (4 $\mu$ g) + A549 whole cell lysate (500 $\mu$ g)

Lane 3: A549 whole cell lysate (20 $\mu$ g)



Immunohistochemistry of paraffin-embedded human kidney tissue using TA386905 at dilution of 1:100



## Western blot

All lanes: DCTN6 antibody at 3.03 µg/ml

Lane 1: Mouse lung tissue

Lane 2: Mouse gonadal tissue

Lane 3: A549 whole cell lysate

Lane 4: U251 whole cell lysate

Lane 5: Mouse thymus tissue

Lane 6: Rat brain tissue

Lane 7: Mouse stomach tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 21 kDa

Observed band size: 21, 18 kDa