

## Product datasheet for **TA386936**


### APOC3 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Apolipoprotein C-III protein (21-99AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">P02656</a>




[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded human liver cancer using TA386936 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using TA386936 at dilution of 1:100



Western blot

All lanes: APOC3 antibody at 2 $\mu$ g/ml + Rat liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 11 kDa

Observed band size: 11 kDa