

Product datasheet for **TA387015M**

GLT8D1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Glycosyltransferase 8 domain-containing protein 1 protein (33-332AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q68CQ7
Background:	Golgi apparatus, membrane, transferase activity, transferring glycosyl groups, polysaccharide biosynthetic process




[View online »](#)


Product images:




Immunohistochemistry of paraffin-embedded human kidney tissue using [TA387015] at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using [TA387015] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human liver tissue using [TA387015] at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using [TA387015] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western blot
All lanes: GLT8D1 antibody at 10 μ g/ml + Mouse liver tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 42, 29 kDa
Observed band size: 42 kDa