

Product datasheet for **TA387021**

IKBIP Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Inhibitor of nuclear factor kappa-B kinase-interacting protein (121-350AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q70UQ0
Background:	Target of p53/TP53 with pro-apoptotic function.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human small intestine tissue using TA387021 at dilution of 1:100



Western blot

All lanes: IKBIPantibody at 4.49μg/ml

Lane 1: Mouse kidney tissue

Lane 2: Mouse lung tissue

Lane 3: Mouse spleen tissue

Lane 4: HL60 whole cell lysate

Lane 5: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 40, 28, 8, 44 kDa

Observed band size: 40, 67, 63 kDa