

Product datasheet for **TA387029**

NR2C2AP Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Nuclear receptor 2C2-associated protein (7-139AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q86WQ0
Background:	May act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes.




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human ovarian cancer using TA387029 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using TA387029 at dilution of 1:100



Western blot

All lanes: NR2C2AP antibody at 4 μ g/ml + HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 16, 20 kDa

Observed band size: 20 kDa