

## Product datasheet for **TA387078M**

### TRIM4 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Reactivity:	Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human E3 ubiquitin-protein ligase TRIM4 protein (1-118AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9C037</a>
Background:	E3 ubiquitin-protein ligase. Mediates 'Lys-63'-linked polyubiquitination of the innate immune receptor DDX58, this linkage doesn't lead to proteasomal degradation but seems to enhance IFN induction.




[View online »](#)


## Product images:



Immunohistochemistry of paraffin-embedded human colon cancer using [TA387078] at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using [TA387078] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



#### Western Blot

Positive WB detected in: MCF-7 whole cell lysate, HeLa whole cell lysate, A549 whole cell lysate, Rat brain tissue, Rat heart tissue

All lanes: TRIM4 antibody at 3.4µg/ml

#### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58, 55, 35 kDa

Observed band size: 55, 62, 35, 90 kDa