

Product datasheet for **TA387080M**

CSRNP2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Cysteine/serine-rich nuclear protein 2 protein (324-543AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9H175
Background:	Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity (By similarity). May play a role in apoptosis.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human testis tissue using [TA387080] at dilution of 1:100



Western blot

All lanes: CSRNP2antibody at 2.34μg/ml

Lane 1: Mouse liver tissue

Lane 2: Mouse heart tissue

Lane 3: Mouse lung tissue

Lane 4: Mouse kidney tissue

Lane 5: HepG2 whole cell lysate

Lane 6: Mouse brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 60 kDa

Observed band size: 60, 48, 24 kDa