

Product datasheet for **TA387178M**

HEATR9 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Protein HEATR9 protein (1-219AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	A2RTY3




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human small intestine tissue using [TA387178] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using [TA387178] at dilution of 1:100



Western blot

All lanes: HEATR9 antibody at 1.2µg/ml + Jurkat whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 66, 27, 61, 45, 40 kDa

Observed band size: 66 kDa