

Product datasheet for **TA387245**


CRIP2 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Cysteine-rich protein 2 protein (1-208AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	P52943
Background:	zinc ion binding



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human pancreatic tissue using TA387245 at dilution of 1:100



Western blot
All lanes: CRIP2 antibody at 4.68 μ g/ml
Lane 1: Rat heart tissue
Lane 2: Mouse lung tissue
Lane 3: MCF-7 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 23, 31 kDa
Observed band size: 23 kDa