

## Product datasheet for **TA387269M**


### OVGP1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Oviduct-specific glycoprotein protein (484-678AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q12889</a>
Background:	Binds to oocyte zona pellucida in vivo. May play a role in the fertilization process and/or early embryonic development.



[View online »](#)

**Product images:**

Western Blot

Positive WB detected in: Mouse heart tissue,  
Mouse kidney tissue, Mouse brain tissue

All lanes: OVGP1 antibody at 3 $\mu$ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 76 kDa

Observed band size: 76 kDa