

Product datasheet for **TA387280**


TNNI3K Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Serine/threonine-protein kinase TNNI3K protein (1-75AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q59H18
Background:	May play a role in cardiac physiology.



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human heart tissue using TA387280 at dilution of 1:100



Western blot
All lanes: TNNI3K antibody at 4.77 μ g/ml
Lane 1: HeLa whole cell lysate
Lane 2: Mouse heart tissue
Lane 3: HepG2 whole cell lysate
Lane 4: Jurkat whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 93, 105, 78, 95 kDa
Observed band size: 93 kDa