

Product datasheet for **TA387383**

RBPJL Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Recombining binding protein suppressor of hairless-like protein (1-151AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9UBG7
Background:	Putative transcription factor, which cooperates with EBNA2 to activate transcription.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human pancreatic cancer using TA387383 at dilution of 1:100



Western Blot
Positive WB detected in: HEK293 whole cell lysate
All lanes: RBPJL antibody at 3 μ g/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 57 kDa
Observed band size: 57 kDa