

## Product datasheet for **TA387578M**

### CYP2F1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Cytochrome P450 2F1 protein (100-290AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">P24903</a>
Background:	May be involved in the metabolism of various pneumotoxins including naphthalene. Is able to dealkylate ethoxycoumarin, propoxycoumarin, and pentoxiresorufin but possesses no activity toward ethoxiresorufin and only trace dearylation activity toward benzyloxiresorufin. Bioactivates 3-methylindole (3MI) by dehydrogenation to the putative electrophile 3-methylene-indolenine.




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human lung tissue using [TA387578] at dilution of 1:100



Western blot  
All lanes: Cytochrome P450 2F1 antibody at 8 µg/ml + A549 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 56, 37 kDa  
Observed band size: 56 kDa