

## Product datasheet for **TA387588**


### CA8 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Carbonic anhydrase-related protein (1-290AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">P35219</a>
Background:	Does not have a carbonic anhydrase catalytic activity.




[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded human brain tissue using TA387588 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using TA387588 at dilution of 1:100



Western blot

All lanes: CA8 antibody at 8 $\mu$ g/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa