

Product datasheet for **TA387777M**


ZNF625 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Zinc finger protein 625 protein (1-306AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q96I27
Background:	May be involved in transcriptional regulation.



[View online »](#)

Product images:

Western blot

All lanes: ZNF625 antibody at 4 μ g/ml

Lane 1: Mouse brain tissue

Lane 2: Mouse lung tissue

Lane 3: Mouse ovarian tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35, 43 kDa

Observed band size: 35 kDa