

Product datasheet for **TA387798**


ADAM33 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 33 protein (140-420AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9BZ11




[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human lung tissue using TA387798 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using TA387798 at dilution of 1:100



Western blot

All lanes: ADAM33 antibody at 4 μ g/ml + NIH/3T3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 88, 85, 33 kDa

Observed band size: 88 kDa