

## Product datasheet for **TA387812M**

### **DENND2D Rabbit Polyclonal Antibody**


#### **Product data:**

|                       |                                                                                                                                                                                                |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                                                                                                                             |
| Applications:         | IHC, WB                                                                                                                                                                                        |
| Recommended Dilution: | Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200                                                                                                                                          |
| Reactivity:           | Human                                                                                                                                                                                          |
| Host:                 | Rabbit                                                                                                                                                                                         |
| Isotype:              | IgG                                                                                                                                                                                            |
| Clonality:            | Polyclonal                                                                                                                                                                                     |
| Immunogen:            | Recombinant Human DENN domain-containing protein 2D protein (1-233AA)                                                                                                                          |
| Formulation:          | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                                                                               |
| Concentration:        | lot specific                                                                                                                                                                                   |
| Purification:         | >95%, Protein G purified                                                                                                                                                                       |
| Conjugation:          | Unconjugated                                                                                                                                                                                   |
| Storage:              | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                  |
| Stability:            | 1 year from dispatch.                                                                                                                                                                          |
| Database Link:        | <a href="#">Q9H6A0</a>                                                                                                                                                                         |
| Background:           | Guanine nucleotide exchange factor (GEF) which may activate RAB9A and RAB9B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form. |




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human small intestine tissue using [TA387812] at dilution of 1:100



Western blot  
All lanes: DENND2D antibody at 4µg/ml + MCF-7 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 54 kDa  
Observed band size: 54 kDa