

## Product datasheet for **TA387839M**

### HAO2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Hydroxyacid oxidase 2 protein (2-200AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9NYQ3</a>
Background:	Catalyzes the oxidation of L-alpha-hydroxy acids as well as, more slowly, that of L-alpha-amino acids.




[View online »](#)


**Product images:**



Immunohistochemistry of paraffin-embedded human kidney tissue using [TA387839] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using [TA387839] at dilution of 1:100



## Western Blot

Positive WB detected in: MCF-7 whole cell lysate, Rat kidney tissue, Mouse kidney tissue

All lanes: HAO2 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39, 41 kDa

Observed band size: 39 kDa