

Product datasheet for **TA387885**


CA12 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Carbonic anhydrase 12 protein (25-301AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	O43570
Background:	Reversible hydration of carbon dioxide.




[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human kidney tissue using TA387885 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using TA387885 at dilution of 1:100



Western blot
All lanes: CA12 antibody at 6 μ g/ml + Mouse
kidney tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 40, 39 kDa
Observed band size: 40, 48 kDa