

## Product datasheet for **TA387889**

### PRR11 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:200-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Proline-rich protein 11 protein (1-360AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q96HE9</a>
Background:	Plays a critical role in cell cycle progression.




[View online »](#)


## Product images:




Immunofluorescent analysis of PC-3 cells using TA387889 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot  
Positive WB detected in: Mouse heart tissue  
All lanes: PRR11 antibody at 3.2µg/ml  
Secondary  
Goat polyclonal to rabbit IgG at 1/50000 dilution  
Predicted band size: 41 kDa  
Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using TA387889 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using TA387889 at dilution of 1:100