


## Product datasheet for **TA387940**

### OXNAD1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Oxidoreductase NAD-binding domain-containing protein 1 protein (18-312AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q96HP4</a>

[View online »](#)

**Product images:****Western blot**

All lanes: OXNAD1 antibody at 4 $\mu$ g/ml

Lane 1: Mouse large intestine tissue

Lane 2: Mouse lung tissue

Lane 3: K562 whole cell lysate

Lane 4: U251 whole cell lysate

Lane 5: A549 whole cell lysate

**Secondary**

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa