

Product datasheet for **TA387990M**

CEACAM4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 4 protein (36-155AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q75871
Background:	Granulocyte orphan receptor that acts as an trigger efficient phagocytosis of attached particles.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human colon cancer using [TA387990] at dilution of 1:100



Western blot
All lanes: CEACAM4 antibody at 1.81µg/ml
Lane 1: Rat heart tissue
Lane 2: Mouse brain tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 26 kDa
Observed band size: 26 kDa