


Product datasheet for **TA388032**

ZNF20 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000–1:5000, IHC:1:20–1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Zinc finger protein 20 protein (1–200AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at –20°C or –80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	P17024
Background:	May be involved in transcriptional regulation.



Product images:

Immunohistochemistry of paraffin-embedded human testis tissue using TA388032 at dilution of 1:100



Western blot

All lanes: ZNF20 antibody at 2.53 μ g/ml + MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 62 kDa

Observed band size: 62, 40 kDa