

## Product datasheet for **TA388062M**

### **BLOC1S1 Rabbit Polyclonal Antibody**


#### **Product data:**

|                       |                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                                |
| Applications:         | IF, IHC, WB                                                                                       |
| Recommended Dilution: | Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200                              |
| Reactivity:           | Human                                                                                             |
| Host:                 | Rabbit                                                                                            |
| Isotype:              | IgG                                                                                               |
| Clonality:            | Polyclonal                                                                                        |
| Immunogen:            | Recombinant Human Biogenesis of lysosome-related organelles complex 1 subunit 1 protein (1-153AA) |
| Formulation:          | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                  |
| Concentration:        | lot specific                                                                                      |
| Purification:         | >95%, Protein G purified                                                                          |
| Conjugation:          | Unconjugated                                                                                      |
| Storage:              | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| Stability:            | 1 year from dispatch.                                                                             |
| Database Link:        | <a href="#">P78537</a>                                                                            |




[View online »](#)


## Product images:



Immunofluorescent analysis of HepG2 cells using [TA388062] at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human small intestine tissue using [TA388062] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using [TA388062] at dilution of 1:100



Western blot  
All lanes: BLOC1S1 antibody at 2 $\mu$ g/ml  
Lane 1: EC109 whole cell lysate  
Lane 2: 293T whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/15000 dilution  
Predicted band size: 18, 15 kDa  
Observed band size: 75 kDa