

Product datasheet for **TA388070M**

MRPL49 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human 39S ribosomal protein L49, mitochondrial protein (1-166AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q13405




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human kidney tissue using [TA388070] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using [TA388070] at dilution of 1:100



Western blot

All lanes: MRPL49 antibody at 0.6µg/ml

Lane 1: Colo320 whole cell lysate

Lane 2: Mouse liver tissue

Lane 3: Mouse kidney tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 20 kDa

Observed band size: 20 kDa