

Product datasheet for **TA388087M**


RNF212 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Probable E3 SUMO-protein ligase RNF212 protein (1-297AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q495C1
Background:	SUMO E3 ligase that acts as a regulator of crossing-over during meiosis: required to couple chromosome synapsis to the formation of crossover-specific recombination complexes. Localizes to recombination sites and stabilizes meiosis-specific recombination factors, such as MutS-gamma complex proteins (MSH4 and MSH5) and TEX11. May mediate sumoylation of target proteins MSH4 and/or MSH5, leading to enhance their binding to recombination sites. Acts as a limiting factor for crossover designation and/or reinforcement and plays an antagonist role with CCNB1IP1/HEI10 in the regulation of meiotic recombination (By similarity).



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human adrenal gland tissue using [TA388087] at dilution of 1:100



Western blot
All lanes: RNF212 antibody at 14 μ g/ml + HepG2 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 34, 15, 22, 32, 27, 31 kDa
Observed band size: 34 kDa