

Product datasheet for **TA388105**

DMKN Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Dermokine protein (1-137AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q6E0U4
Background:	May act as a soluble regulator of keratinocyte differentiation.



[View online »](#)

Product images:



Western blot

All lanes: DMKN antibody at 6 μ g/ml

Lane 1: Mouse kidney tissue

Lane 2: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 16, 37, 46, 39, 42, 36, 21,
19, 17, 15, 10, 47 kDa

Observed band size: 42 kDa