

Product datasheet for **TA388122**


RAD9B Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Cell cycle checkpoint control protein RAD9B protein (274-408AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q6WBX8
Background:	checkpoint clamp complex, nucleoplasm, 3'-5' exonuclease activity, cellular response to ionizing radiation, DNA repair, DNA replication, DNA replication checkpoint, intra-S DNA damage checkpoint, nucleic acid phosphodiester bond hydrolysis, regulation of signal transduction by p53 class mediator



[View online »](#)

Product images:


Western Blot

Positive WB detected in: Mouse brain tissue

All lanes: RAD9B antibody at 6 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 48, 47, 39, 29, 30 kDa

Observed band size: 48 kDa