

## Product datasheet for **TA388197M**

### ISOC1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Isochorismatase domain-containing protein 1 protein (155-298AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q96CN7</a>
Background:	cytoplasm, extracellular exosome, peroxisome




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human breast cancer using [TA388197] at dilution of 1:100



Western Blot  
Positive WB detected in: SH-SY5Y whole cell lysate  
All lanes: ISOC1 antibody at 3.5µg/ml  
Secondary  
Goat polyclonal to rabbit IgG at 1/50000 dilution  
Predicted band size: 33 kDa  
Observed band size: 33 kDa