

Product datasheet for **TA388198**

RMDN1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Regulator of microtubule dynamics protein 1 protein (13-314AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q96DB5




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using TA388198 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using TA388198 at dilution of 1:100



Western blot

All lanes: RMDN1 antibody at 0.6 μ g/ml

Lane 1: K562 whole cell lysate

Lane 2: HepG2 whole cell lysate

Lane 3: MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 36, 33, 31 kDa

Observed band size: 31 kDa