


Product datasheet for **TA388218**

TSGA13 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500–1:2000, IHC:1:20–1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Testis-specific gene 13 protein (1–275AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at –20°C or –80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q96PP4
Background:	positive regulation of protein targeting to mitochondrion




Product images:

Immunohistochemistry of paraffin-embedded human testis tissue using TA388218 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using TA388218 at dilution of 1:100



Western blot

All lanes: TSGA13 antibody at 4 μ g/ml + A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 32 kDa

Observed band size: 32 kDa