

Product datasheet for **TA388232M**

UTP23 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human rRNA-processing protein UTP23 homolog protein (1-249AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9BRU9
Background:	Involved in rRNA-processing and ribosome biogenesis.




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human liver cancer using [TA388232] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using [TA388232] at dilution of 1:100



Western blot
 All lanes: UTP23 antibody at 7µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 29, 17 kDa
 Observed band size: 29 kDa