

## Product datasheet for **TA388235**

### **NOL10 Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Nucleolar protein 10 protein (1-210AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9BSC4</a>
Background:	nucleolus, poly(A) RNA binding




[View online »](#)


## Product images:



Immunohistochemistry of paraffin-embedded human pancreatic tissue using TA388235 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using TA388235 at dilution of 1:100



Western blot

All lanes: NOL10 antibody at 8µg/ml + K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 81, 75, 22, 78 kDa

Observed band size: 81 kDa