

Product datasheet for **TA388250**


ROGDI Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Protein rogdi homolog protein (14-287AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q9GZN7
Background:	May act as a positive regulator of cell proliferation.




[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human kidney tissue using TA388250 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using TA388250 at dilution of 1:100



Western blot

All lanes: ROGDI antibody at 1.2 μ g/ml

Lane 1: HeLa whole cell lysate

Lane 2: K562 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa