

## Product datasheet for **TA388280**

### **DNAJC17 Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Dnaj homolog subfamily C member 17 protein (1-304AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9NVM6</a>




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human colon cancer using TA388280 at dilution of 1:100



Western blot  
All lanes: DNAJC17 antibody at 16µg/ml  
Lane 1: A431 whole cell lysate  
Lane 2: Hela whole cell lysate  
Lane 3: Raji whole cell lysate  
Lane 4: HepG2 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 35 kDa  
Observed band size: 35 kDa