

Product datasheet for **TA388435M**

GGTLC2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Glutathione hydrolase light chain 2 protein (1-218AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q14390
Background:	anchored component of external side of plasma membrane, extracellular exosome, gamma-glutamyltransferase activity, glutathione metabolic process, leukotriene biosynthetic process



[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human liver cancer using [TA388435] at dilution of 1:100



Western blot
All lanes: GGTLC2 antibody at 2µg/ml + Mouse kidney tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 24 kDa
Observed band size: 24 kDa