

Product datasheet for **TA388477M**

DHRS13 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Dehydrogenase/reductase SDR family member 13 protein (1-327AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q6UX07
Background:	Putative oxidoreductase.




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human brain tissue using [TA388477] at dilution of 1:100



Western Blot
Positive WB detected in: HEK293 whole cell lysate
All lanes: DHRS13 antibody at 4µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 41, 36, 33 kDa
Observed band size: 41 kDa