

## Product datasheet for **TA388502**

### ODF3L1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Outer dense fiber protein 3-like protein 1 protein (1-274AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q8IXM7</a>




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human gastric cancer using TA388502 at dilution of 1:100



Western blot  
All lanes: ODF3L1 antibody at 3µg/ml + HT29 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 32 kDa  
Observed band size: 32 kDa