

Product datasheet for **TA388533**

SUMF2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IP, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IP:1:1000-1:2000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Sulfatase-modifying factor 2 protein (26-301AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Antigen Affinity purified & Affinity purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q8NBJ7
Background:	Lacks formyl-glycine generating activity and is unable to convert newly synthesized inactive sulfatases to their active form. Inhibits the activation of sulfatases by SUMF1.




[View online »](#)

Product images:



Immunoprecipitating SUMF2 in HeLa whole cell lysate
 Lane 1: Rabbit control IgG instead of TA388533 in HeLa whole cell lysate. For western blotting, a HRP-conjugated Protein G antibody was used as the secondary antibody (1/50000)
 Lane 2: TA388533 (5µg) + HeLa whole cell lysate (500µg)
 Lane 3: HeLa whole cell lysate (10µg)



Western Blot
 Positive WB detected in: A549 whole cell lysate, HepG2 whole cell lysate, Jurkat whole cell lysate, K562 whole cell lysate, U87 whole cell lysate, HeLa whole cell lysate
 All lanes: SUMF2 antibody at 1:2000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 34, 25, 38, 24, 35 kDa
 Observed band size: 34 kDa