

Product datasheet for **TA388553M**

RUFY2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human RUN and FYVE domain-containing protein 2 protein (396-655AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<u>Q8WXA3</u>




[View online »](#)


Product images:



Immunohistochemistry of paraffin-embedded human testis tissue using [TA388553] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using [TA388553] at dilution of 1:100



Western blot

All lanes: RUFY2 antibody at 2.4µg/ml + PC-3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 76, 71, 75, 47, 44 kDa

Observed band size: 76, 55, 49 kDa