

Product datasheet for **TA388592**


PTPN18 Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Reactivity:	Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Tyrosine-protein phosphatase non-receptor type 18 protein (1-210AA)
Formulation:	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q99952
Background:	Differentially dephosphorylate autophosphorylated tyrosine kinases which are known to be overexpressed in tumor tissues.



[View online »](#)

Product images:

Immunohistochemistry of paraffin-embedded human brain tissue using TA388592 at dilution of 1:100



Western blot
All lanes: PTPN18 antibody IgG at 2.98 μ g/ml + Rat heart tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 51, 39 kDa
Observed band size: 51 kDa