

## Product datasheet for **TA388599M**

### FGFBP2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Fibroblast growth factor-binding protein 2 protein (20-223AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9BYJ0</a>




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human testis tissue using [TA388599] at dilution of 1:100



Western blot  
All lanes: FGFBP2 antibody at 2µg/ml + Human serum  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 25 kDa  
Observed band size: 25 kDa