

## Product datasheet for **TA388624M**

### NRBP2 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, IP
Recommended Dilution:	Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200, IP:1:200-1:2000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human Nuclear receptor-binding protein 2 protein (12-258AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9NSY0</a>
Background:	May regulate apoptosis of neural progenitor cells during their differentiation.




[View online »](#)


## Product images:




Immunofluorescence staining of HeLa cells with [TA388624] at 1:133, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunoprecipitating NRBP2 in A549 whole cell lysate  
 Lane 1: Rabbit control IgG (1µg) instead of [TA388624] in A549 whole cell lysate. For western blotting, a HRP-conjugated Protein G antibody was used as the secondary antibody (1/2000)  
 Lane 2: [TA388624] (6µg) + A549 whole cell lysate (500µg)  
 Lane 3: A549 whole cell lysate (10µg)



Immunohistochemistry of paraffin-embedded human endometrial cancer using [TA388624] at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using [TA388624] at dilution of 1:100