

## Product datasheet for **TA388627M**

### OSER1 Rabbit Polyclonal Antibody


#### Product data:

|                       |                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                    |
| Applications:         | IHC, WB                                                                               |
| Recommended Dilution: | Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200                                |
| Reactivity:           | Mouse, Human                                                                          |
| Host:                 | Rabbit                                                                                |
| Isotype:              | IgG                                                                                   |
| Clonality:            | Polyclonal                                                                            |
| Immunogen:            | Recombinant Human Oxidative stress-responsive serine-rich protein 1 protein (1-292AA) |
| Formulation:          | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4      |
| Concentration:        | lot specific                                                                          |
| Purification:         | >95%, Protein G purified                                                              |
| Conjugation:          | Unconjugated                                                                          |
| Storage:              | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| Stability:            | 1 year from dispatch.                                                                 |
| Database Link:        | <a href="#">Q9NX31</a>                                                                |




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human small intestine tissue using [TA388627] at dilution of 1:100



Western blot  
All lanes: OSER1 antibody at 3µg/ml + Mouse liver tissue  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 32 kDa  
Observed band size: 32 kDa