

## Product datasheet for **TA388631**


### ABI3 Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Mouse, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human ABI gene family member 3 protein (1-366AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q9P2A4</a>
Background:	May inhibit tumor metastasis (By similarity). In vitro, reduces cell motility.



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded human tonsil tissue using TA388631 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using TA388631 at dilution of 1:100



Western blot

All lanes: ABI3 antibody at 5 $\mu$ g/ml + Mouse heart tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 40, 39 kDa

Observed band size: 40 kDa