

## Product datasheet for **TA388658**

### **RUNDC3A Rabbit Polyclonal Antibody**


#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human RUN domain-containing protein 3A protein (1-405AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	Antigen Affinity Purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	<a href="#">Q59EK9</a>
Background:	May act as an effector of RAP2A in neuronal cells.




[View online »](#)

## Product images:



Immunohistochemistry of paraffin-embedded human testis tissue using TA388658 at dilution of 1:100



Western blot  
All lanes: RUNDC3A antibody at 1.5µg/ml + MCF-7 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 50, 46, 42 kDa  
Observed band size: 50 kDa