

## Product datasheet for **TA388753**


### **FAM53B Rabbit Polyclonal Antibody**

#### **Product data:**


<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	IHC, WB
<b>Recommended Dilution:</b>	Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
<b>Reactivity:</b>	Human
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	Recombinant Human Protein FAM53B protein (1-305AA)
<b>Formulation:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Concentration:</b>	lot specific
<b>Purification:</b>	>95%, Protein G purified
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Stability:</b>	1 year from dispatch.
<b>Database Link:</b>	<a href="#">Q14153</a>
<b>Background:</b>	Acts as a regulator of Wnt signaling pathway by regulating beta-catenin (CTNNB1) nuclear localization.



[View online »](#)

**Product images:**

Immunohistochemistry of paraffin-embedded human breast cancer using TA388753 at dilution of 1:100



Western blot  
All lanes: FAM53B antibody at 2 $\mu$ g/ml  
Lane 1: MCF-7 whole cell lysate  
Lane 2: A549 whole cell lysate  
Lane 3: HCT116 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 46, 34 kDa  
Observed band size: 46 kDa