

Product datasheet for **TA388783**

OBSL1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human OBSL1 protein (1-61AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q6P710




[View online »](#)

Product images:



Immunohistochemistry of paraffin-embedded human placenta tissue using TA388783 at dilution of 1:100



Western Blot
Positive WB detected in: HepG2 whole cell lysate, MCF-7 whole cell lysate, A549 whole cell lysate,
All lanes: OBSL1 antibody at 2.5µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 207 kDa
Observed band size: 207 kDa