

Product datasheet for **TA388791**


WDR88 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	Recommended dilution: WB:1:500–1:2000, IHC:1:20–1:200
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant Human WD repeat-containing protein 88 protein (1–250AA)
Formulation:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Concentration:	lot specific
Purification:	>95%, Protein G purified
Conjugation:	Unconjugated
Storage:	Upon receipt, store at –20°C or –80°C. Avoid repeated freeze.
Stability:	1 year from dispatch.
Database Link:	Q6ZMY6



Product images:



Immunohistochemistry of paraffin-embedded human placenta tissue using TA388791 at dilution of 1:100



Western blot

All lanes: WDR88 antibody at 6 μ g/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 53, 48 kDa

Observed band size: 53 kDa